
Product Datasheet

PTK2 (phospho-Tyr576+Tyr577) antibody (orb5220)

Description

PTK2 (phospho-Tyr576+Tyr577)

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FAK around the phosphorylation site of Tyr576/577 (ST(p-Y)(p-Y)KA)
Target	phospho-PTK2 (Tyr576+Tyr577)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	119 kDa
Uniprot ID	Q05397
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.