

Product Datasheet

TRIM33 antibody (Biotin) (orb52229)

| | |
|----------------------------|--|
| Description | TRIM33 antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Biotin |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human E3 ubiquitin-protein ligase TRIM33 protein (1-273AA) |
| Target | TRIM33 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9UPN9 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |