


Product Datasheet

SIRT5 antibody (orb522905)

| | |
|----------------------------|---|
| Description | SIRT5 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Synthetic peptide of Human SIRT5 |
| Target | SIRT5 |
| Preservatives | -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9NXA8 |
| Dilution Range | ELISA:1:1000-1:5000, WB:1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Mouse heart ti...