


Product Datasheet

CCZ1 antibody (orb52410)

| | |
|----------------------------|---|
| Description | CCZ1 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant Human Vacuolar fusion protein CCZ1 homolog protein (201-482AA) |
| Target | CCZ1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P86791 |
| Dilution Range | WB:1:500-5000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Rat liver(lane ...