

Product Datasheet

Exendin 4 antibody (orb524139)

Description	Exendin 4 antibody
Species/Host	Rabbit
Reactivity	Reptile
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic peptide of Heloderma suspectum Exendin 4
Target	Exendin-4
Preservatives	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P26349
Expiration Date	12 months from date of receipt.