



Product Datasheet

FLT3 (phospho-Tyr842) antibody (orb5253)


Description	FLT3 (phospho-Tyr842) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FLT3 around the phosphorylation site of Tyr842 (SN(p-Y)VV)
Target	Phospho-FLT3 (Tyr842)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	113 kDa
Uniprot ID	P36888
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat lung...



Immunohistochemical staining of rat brai...



Immunohistochemical staining of Rat lung...