

Product Datasheet

ML RR-S2 CDA (orb525605)

Description	ML RR-S2 CDA is a synthetic cyclic dinucleotide (CDN) that contains non-canonical 2'5'-phosphodiester bonds and is an activator of stimulator of interferon genes (STING).
Note	For research use only
Formula	$C_{20}H_{24}N_{10}O_{10}P_2S_2$
MW	690.54
CAS Number	[1638241-89-0]
SMILES	<chem>NC1=NC=NC2=C1N=CN2[C@@H]3O[C@@H]4CO[P@](S)(=O)O[C@@H]5[C@H](O)[C@@H](CO[P@](S)(=O)O[C@H]4[C@H]3O)O[C@H]5N6C=NC7=C6N=CN=C7N</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)