

Product Datasheet

NGI-1 (ML414) (orb525636)

Description	NGI-1 (ML414) is a potent, cell-permeable oligosaccharyltransferase (OST) inhibitor with IC50 of 1.1...
Note	For research use only
Formula	$C_{17}H_{22}N_4O_3S_2$
MW	394.51
CAS Number	[790702-57-7]
SMILES	<chem>O=C(NC1=NC=C(C)S1)C2=CC(S(=O)(N(C)C)=O)=CC=C2N3CCCC3</chem>
Expiration Date	12 months from date of receipt.