

Product Datasheet

Flag Tag antibody (orb526284)

Description	Flag Tag antibody
Species/Host	Mouse
Reactivity	All
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Flag Tag
Target	Flag Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB=1:5000-20000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.