

Product Datasheet

ERK1 (phospho-Thr202) antibody (orb526443)

Description	Rabbit polyclonal antibody to ERK1 (phospho-Thr202) conjugated to HRP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse ERK1/2 around the phosphorylation site of Thr202/Tyr204 FL(p-T)E(p-Y)VA
Target	MAPK3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q63844
Dilution Range	ELISA: 1:100-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.