

Product Datasheet

Raf1 (phospho-Ser43) antibody (orb526683)

Description	Raf1 (phospho-Ser43) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human RAF1 around the phosphorylation site of Ser43 (RA(p-S)DD)
Target	Phospho-RAF1 (Ser43)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	71 kDa
Uniprot ID	P04049
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.