


## Product Datasheet

# Human CA1/Carbonic anhydrase 1 ELISA Kit (orb526784)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human CA1/Carbonic anhydrase 1 ELISA Kit   |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 312 pg/ml - 20,000 pg/ml   |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Expression system for standard: E.coli;<br>Immunogen sequence: A2-F261   |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months.<br>Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Dilution ratio of 1:100, concentration in serum and plasma is around 600 ng/ml   |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum, plasma (heparin, EDTA, citrate) and urine.   |
| <b>Uniprot ID</b>          | <b>P00915</b>  |
| <b>Sensitivity</b>         | 50 pg/ml   |
| <b>Expiration Date</b>     | Please enquire.  |



Human CA1  
PicoKine  
ELISA  
Kitstandard  
cur...