


## Product Datasheet

### Rat IL6ST ELISA Kit (orb526796)

<b>Description</b>	Rat IL6ST ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.78 ng/ml - 50 ng/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: Q23-E618
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:10, concentration in serum and plasma is around 40 ng/ml
<b>Sample Types</b>	cell culture supernatants, serum and plasma (heparin, EDTA, citrate ).
<b>Uniprot ID</b>	<b>P40190</b>
<b>Sensitivity</b>	0.05 ng/ml
<b>Expiration Date</b>	Please enquire.



Rat IL6ST  
PicoKine  
ELISA Kit  
standard  
CU...