


Product Datasheet

Mouse TNFRSF12A/TWEAKR ELISA Kit (orb526801)

Description	Mouse TNFRSF12A/TWEAKR ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E28-P80
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 2 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q9CR75
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse
TNFRSF12A/TWEAKR
PicoKine ELISA Ki...