


Product Datasheet

Mouse LYPD3 ELISA Kit (orb526826)

Description	Mouse LYPD3 ELISA Kit
Reactivity	Mouse
Range	93.6 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L33-S343
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 200-2000 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q91YK8
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse
LYPD3
PicoKine
ELISA Kit
standard ...