


Product Datasheet

Human MANF ELISA Kit (orb526831)

Description	Human MANF ELISA Kit
Reactivity	Human
Range	93.6 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L25-L182
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:5, concentration in serum and plasma is around 8 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P55145
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Human
MANF
PicoKine
ELISA Kit
standard c...