


Product Datasheet

Human Serpin A6 ELISA Kit (orb526880)

Description	Human Serpin A6 ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M23-V405
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is 0.5-1 ug/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA) and urine.
Uniprot ID	P08185
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Human
Serpin A6
PicoKine
ELISA
Kitstanda...