


Product Datasheet

Human FLRT2 ELISA Kit (orb526886)

Description	Human FLRT2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C36-T660
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 50 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	O43155
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
FLRT2
PicoKine
ELISA
Kitstandard
C...