


Product Datasheet

Mouse CRELD2 ELISA Kit (orb526898)

Description	Mouse CRELD2 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S23-L350
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:4, concentration in serum and plasma is around 6 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q9CYA0
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse
CRELD2
PicoKine
ELISA Kit
Standard...