



Product Datasheet

Aromatase/CYP19A1 Antibody (orb526990)


Description	Aromatase/CYP19A1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC, IHC-Fr, WB
Immunogen	E. coli-derived human CYP19A1 recombinant protein (Position: Y241-H503).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	P11511
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of U2OS cells us...



WB analysis of Aromatase using anti-Arom...



IHC analysis of Aromatase using anti-Aro...