

Product Datasheet

Hyaluronan synthase 1/HAS1 Antibody (orb526995)

Description

Hyaluronan synthase 1/HAS1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence of human HAS1 (NRAEDLYMVDMFREVFADEDPATYVWDGNYHQWPEA)

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

70 kDa

Uniprot ID


[Q92839](#)

Dilution Range


Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunohistochemistry (Frozen Section), 0.5-1µg/ml Immunocytochemistry/Immunofluorescence, 2µg/ml Flow Cytometry, 1-3µg/1x10⁶ cells

Expiration Date


12 months from date of receipt.



Flow
Cytometry
analysis of U-
87MG cells ...



WB:Lane
1:SHG-44
cell;2:THP-1
cell;3:rat...



IF analysis of
HAS1 using
anti-HAS1
anti...