


Product Datasheet


uPA Receptor/PLAUR Antibody (orb527001)

Description	uPA Receptor/PLAUR Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human uPA Receptor/PLAUR (TKSGCNHPDLDVQYRS).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q03405
Expiration Date	12 months from date of receipt.




Flow Cytometry analysis of H-PBMC cells

...



Flow Cytometry analysis of H-PBMC cells

...



WB analysis of uPA Receptor using anti-uPA antibody