



Product Datasheet

PKC beta 1/PRKCB Antibody (orb527003)


Description	PKC beta 1/PRKCB Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human PKC beta 1/PRKCB (NFDKEFTRQPVELTPTDKLFIMNLDQNE).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	77 kDa
Uniprot ID	P05771
Expiration Date	12 months from date of receipt.



WB analysis of PKC beta 1 using anti-PKC...



WB analysis of PKC beta 1 using anti-PKC...



IHC analysis of PKC beta 1 using anti-PK...