



Product Datasheet

WWOX Antibody (monoclonal, 3D10) (orb527044)


Description	WWOX Antibody (monoclonal, 3D10)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E. coli-derived human WWOX recombinant protein (Position: M1-D245).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	47 kDa
Uniprot ID	Q9NZC7
Sensitivity	10pg/ml
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of U2OS cells us...



WB:1:human placenta tissue;2:T-47D;3:Hep...



IF analysis of WWOX using anti-WWOX anti...