
Product Datasheet

**GM130 GOLGA2 Antibody (monoclonal, 4G3)
(orb527062)**

Description

GM130 GOLGA2 Antibody (monoclonal, 4G3)


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications

Flow Cytometry analysis of A431 cells us...

Immunogen

E.coli-derived human GM130 recombinant protein (Position: E796-E913).

Form/Appearance


Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.



WB analysis of GM130 using anti-GM130 an...

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500µg/ml.

Isotype

Mouse IgG2b

Clonality


Monoclonal

MW

140 kDa

Uniprot ID
[Q08379](https://www.uniprot.org/entry/Q08379)
Expiration Date

12 months from date of receipt.



IF analysis of GM130 using anti-GM130 an...