

## Product Datasheet

### Chenodeoxycholic acid-BSA conjugate (orb527307)

<b>Description</b>	Chenodeoxycholic acid-BSA conjugate
<b>Tested Applications</b>	ELISA, IHC
<b>Preservatives</b>	0.01mol/L, pH 7.4 PBS
<b>Form/Appearance</b>	Dry powder
<b>Storage</b>	Short-term store at 4°C. Long-term store at -23°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Application notes</b>	Other Small Molecule Antigen
<b>Source</b>	Bovine serum albumin (BSA)
<b>CAS Number</b>	474-25-9
<b>Hazard Information</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.