

## Product Datasheet

### Retinoic acid- HSA conjugate (orb527386)

<b>Description</b>	Retinoic acid- HSA conjugate
<b>Tested Applications</b>	ELISA, ICA
<b>Preservatives</b>	0.01mol/L, pH 7.4 PBS
<b>Form/Appearance</b>	Dry powder
<b>Storage</b>	Short-term store at 4°C. Long-term store at -20°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Application notes</b>	Other Small Molecule Antigen
<b>Source</b>	Human serum albumin (HSA)
<b>CAS Number</b>	302-79-4
<b>Hazard Information</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.