

## Product Datasheet

### Mouse Immunoglobulin G1 protein (orb527408)

<b>Description</b>	Mouse Immunoglobulin G1 protein
<b>Tested Applications</b>	ELISA, SDS-PAGE, WB
<b>Target</b>	Mouse Immunoglobulin G1
<b>Preservatives</b>	PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	> 95%(SDS-PAGE)
<b>Source</b>	Purified from pooled normal mouse serum
<b>Hazard Information</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.