

Product Datasheet

Mouse Immunoglobulin G2A protein (orb527409)

Description	Mouse Immunoglobulin G2A protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Mouse Immunoglobulin G2a
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from pooled normal mouse serum
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.