

Product Datasheet

Goat IgG Fc Fragment protein (orb527412)

| | |
|----------------------------|---|
| Description | Goat IgG Fc Fragment protein |
| Tested Applications | ELISA, SDS-PAGE, WB |
| Target | Goat IgG Fc fragment |
| Preservatives | PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Storage | Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Purity | > 95%(SDS-PAGE) |
| Source | Goat serum IgG digested with papain |
| Hazard Information | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)