

Product Datasheet

Rabbit IgG Fc Fragment protein (orb527415)

Description	Rabbit IgG Fc Fragment protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Rabbit IgG Fc fragment
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Rabbit serum IgG
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.