

Product Datasheet

Rat IgG Fc Fragment protein (orb527416)

Description	Rat IgG Fc Fragment protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Rat IgG Fc fragment
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Rat serum IgG digested with papain
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.