

Product Datasheet

Rat IgG Fab Fragment protein (orb527421)

| | |
|----------------------------|---|
| Description | Rat IgG Fab Fragment protein |
| Tested Applications | ELISA, SDS-PAGE, WB |
| Target | Rat IgG Fab fragment |
| Preservatives | PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Storage | Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Purity | > 95%(SDS-PAGE) |
| Source | Rat serum IgG digested with papain |
| Hazard Information | For research use only |
| Expiration Date | 6 months from date of receipt. |