

Product Datasheet

Human Neutrophil Gelatinase-Associated Lipocalin Protein (orb527431)

Description	Human Neutrophil Gelatinase-Associated Lipocalin Protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	human Neutrophil gelatinase-associated lipocalin protein
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from human urine
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.