

Product Datasheet

Canine Immunoglobulin G protein (orb527434)

Description	Canine Immunoglobulin G protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Canine Immunoglobulin G
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from pooled normal Canine serum
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.