

## Product Datasheet

### Guinea Pig Immunoglobulin G protein (orb527437)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Guinea Pig Immunoglobulin G protein  |
| <b>Tested Applications</b> | ELISA, SDS-PAGE, WB  |
| <b>Target</b>              | Guinea pig Immunoglobulin G  |
| <b>Preservatives</b>       | PBS, pH 7.4  |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Aliquot and store at -20°C or -80°C.<br>Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Purity</b>              | > 95%(SDS-PAGE)  |
| <b>Source</b>              | Purified from pooled normal guinea pig serum                               |
| <b>Hazard Information</b>  | For research use only  |
| <b>Expiration Date</b>     | 6 months from date of receipt.   |