

Product Datasheet

Rabbit Immunoglobulin G protein (orb527443)

| | |
|----------------------------|--|
| Description | Rabbit Immunoglobulin G protein |
| Tested Applications | ELISA, SDS-PAGE, WB |
| Target | Rabbit Immunoglobulin G |
| Preservatives | PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Storage | Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Purity | > 95%(SDS-PAGE) |
| Source | Purified from Rabbit plasma |
| Hazard Information | For research use only |
| Expiration Date | 6 months from date of receipt. |