

Product Datasheet

Wheat Gliadin protein (orb527468)

Description	Wheat Gliadin protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Wheat Gliadin
Preservatives	70% Ethanol
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from wheat seed
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.