

Product Datasheet

Immobilized Papain protein (orb527473)

Description	Immobilized Papain protein
Tested Applications	Functional Studies
Target	Immobilized Papain
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from papaya latex
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.