

Product Datasheet

SLK (phospho-Tyr530) antibody (orb5277)

Description	SLK (phospho-Tyr530) antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FYN around the phosphorylation site of Tyr530 (PQ(p-Y)QP)
Target	phospho-SLK (Tyr530)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	P06241
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.