

## Product Datasheet

### PIK3R2 antibody (Biotin) (orb52858)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | PIK3R2 antibody (Biotin)   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Biotin   |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Recombinant Human Phosphatidylinositol 3-kinase regulatory subunit beta protein (75-125AA) |
| <b>Target</b>              | PIK3R2   |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                             |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                              |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/O00459">O00459</a>                                  |
| <b>Dilution Range</b>      | ELISA  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |