

Product Datasheet

POLR1C antibody (FITC) (orb52877)

Description	POLR1C antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human DNA-directed RNA polymerases I and III subunit RPAC1 protein (125-230AA)
Target	POLR1C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15160
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.