

Product Datasheet

GATA6 (phospho-Tyr271) antibody (orb5304)

Description	GATA6 (phospho-Tyr271) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GATA6 around the phosphorylation site of Tyr271 (GL(p-Y)SK)
Target	phospho-GATA6 (Tyr271)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60 kDa
Uniprot ID	Q92908
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.