

Product Datasheet

Probable terminase, ATPase subunit antibody (orb53259)

Description	Probable terminase, ATPase subunit antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Haemophilus phage HP1 Probable terminase, ATPase subunit protein (1-190AA)
Target	Probable terminase, ATPase subunit
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51718
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.