


Product Datasheet

petE antibody (orb53267)

Description	petE antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Prochlorothrix hollandica Plastocyanin protein (35-131AA)
Target	petE
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P50057
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...