

Product Datasheet

Mid Range DNA Ladder (orb532917)

Description	Mid Range DNA Ladder consists of 12 DNA fragments ranging from 100 bp to 3 kb. It is designed to sho...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 8% (w/v) Glycerol, Tartrazine and Xylene Cyanol FF; Color: green
Concentration	130 ng/μl
Storage	store at -20 °C
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.