

## Product Datasheet

### 50 bp DNA Ladder (orb532921)

<b>Description</b>	50 bp DNA Ladder consists of 10 DNA fragments ranging from 50 bp to 500 bp and is supplied in ready-...
<b>Form/Appearance</b>	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol and Orange G; Color: orange
<b>Concentration</b>	105 ng/μl
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Application notes</b>	<b>Recommended Load:</b> 5 μl per lane.
<b>Expiration Date</b>	12 months from date of receipt.