

Product Datasheet

200 bp DNA Ladder (orb532923)

| | |
|--------------------------|---|
| Description | 200 bp DNA Ladder consist of 10 DNA fragments ranging from 200 bp to 2 kb and is supplied in ready-t... |
| Form/Appearance | liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, Bromophenol Blue and Xylene Cyanol FF; Color: blue |
| Concentration | 105 ng/μl |
| Storage | store at -20 °C |
| Note | For research use only |
| Application notes | Recommended Load: 5 μl per lane. |
| Expiration Date | 12 months from date of receipt. |