

Product Datasheet

Fluorescent High Range DNA Ladder (orb532928)

| | |
|--------------------------|--|
| Description | Fluorescent High Range DNA Ladder consists of 12 DNA fragments ranging from 500 bp to 10 kb. It is d... |
| Form/Appearance | liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Green, Bromophenol Blue and Xylene Cyanol FF; Color: blue |
| Concentration | 130 ng/μl |
| Storage | store at -20 °C |
| Note | For research use only |
| Application notes | Recommended Load: 5 μl per lane. |
| Expiration Date | 12 months from date of receipt. |