

## Product Datasheet

### Fluorescent 100 bp DNA Ladder (orb532932)

<b>Description</b>	Fluorescent 100 bp DNA Ladder consist of 10 DNA fragments ranging from 100 bp to 1 kb. The prestaine...
<b>Form/Appearance</b>	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Tartrazine and Xylene Cyanol FF; Color: green
<b>Concentration</b>	105 ng/μl
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Application notes</b>	<b>Recommended Load:</b> 5 μl per lane.
<b>Expiration Date</b>	12 months from date of receipt.